




neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLab® Screwdriver Set 1000 V

€58.40
plus VAT &
Shipping

Product Images



Description

High quality tool set with blades made of chrome vanadium molybdenum steel.

Insulationsprayed directly onto the blade.

The set consists of:

1 each slotted screwdriver: 2.8 x 0.6 mm, 4.0 x 0.8 mm, 5.5 x 1 mm and 6.5 x 1.2 mm.

1 each Phillips screwdriver: size 1 and 2.

1 voltage tester.

Can be used up to 1000 volts.

VDE and GS mark, unit-tested.

Protective insulation according to DIN EN 60900.

Additional Information

No.	2-2007
Manufacturer (Brand)	neoLab
EAN	4058072059330
Transport temperature	Room temperature
Material	Steel chrome plated
Tools properties	in set
Functionality	mechanical
Number of parts in set	7
Width slotted screwdriver	2.8 mm 4.0 mm 5.5 mm 6.5 mm
Size Phillips screwdriver	1 mm 2 mm
Strength slotted screwdriver	0.6 mm 0.8 mm 1.0 mm 1.2 mm
Type Screwdriver / Wrench	Phillips screwdriver Slotted head screwdriver

