




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® BagOpen for BagFilter and BagPage

€64.40
plus VAT &
Shipping

Product Images



Description

BagOpen® automatically opens bags and holds them open during sample insertion or removal. Four self-adhesive rubbers adhere to the outer walls of the bag and hold it open. This avoids the risk of coming into contact with the sample.

- for bags with 400 ml volume
- autoclavable (except rubbers)
- 4 freely positionable rubbers incl.

Additional Information

No.	2-1385
Manufacturer (Brand)	neoLab
EAN	4058072056377
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Material	Polycarbonate (PC)
for type bag	BagFilter 400 ml for BagMixer BagPage 400 ml for BagMixer Sterile bags for Bagmixer 400VW

