



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Lid opener plastic non-slip, 11 cm Ø, 2 pcs./pack

€19.50
plus VAT &
Shipping

Product Images



Description

Soft plastic cone for opening lids on wide-necked bottles, jars, etc.

The non-slip material ensures the right grip even when closing.

11 cm diameter, 4 mm thick.

Temperature resistant from 0 to 50°C.

Additional Information

No.	3-0924
Manufacturer (Brand)	neoLab
EAN	4058072011734
Transport temperature	Room temperature
diameter	110 mm
Color	blue
Material	Silicone
TBST MAX	50 °C
TBST MIN	0 °C
Thickness	4 mm

