




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® cable tie pliers for ties 2.2 to 4.8 mm

€31.60
plus VAT &
Shipping

Product Images



Description

Cable tie locking pliers, infinitely adjustable for universal ties from 2.2 to 4.8 mm wide.

Mechanical parts made of steel, housing made of plastic.

With adjustable clamping force regulation, automatic cutter.

With large clamping lever for easy tightening of cable ties.

Can be used e.g. for 1-7110 ff.

Additional Information

No.	1-7115
Manufacturer (Brand)	neoLab
EAN	4058072045166
Color	red
Material	Steel
Functionality	mechanical
Type pliers	Locking pliers

