




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hecht hemacytometer cover slips 22 x 22 mm, 10 pcs/pack

€13.10 As low as
€12.00
plus VAT &
Shipping

Product Images



Description

Cover glasses for hemacytometers (blood analysis), precision design, with optically flat ground surfaces on both sides. CE design. Made of borosilicate glass, according to DIN ISO 8255.

Refractive index_{ne} = 1.52 ± 0.01.

Thickness: 4

Thickness: 0.4 mm

PU: 10 pcs./pack

Additional Information

No.	WQ-0307
Manufacturer (Brand)	Hecht
old neoLab article no.	C-1006
EAN	4058072443436
Transport temperature	Room temperature
Material	Borosilicate glass
Wide	22 mm
Form	angular
Length	22 mm
Application cover glass	for hemacytometer
Thickness	0.4 mm
Dimensions	22 x 22 mm
Strength/Class	4

