




neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLab® slide holder for staining box

€24.10
plus VAT &
Shipping

Product Images



Description

The staining rack with 12 individual places for slides (75 x 25 or 76 x 26 mm).

The labeling fields are not stained as well, so that removal is possible without tweezers. The lid of the staining rack is flush with the staining box and has a handle. When the lid is placed over the slides, they are optimally protected. The lid can be moved sideways for removal.

Temperature resistant from -170°C to +121°C, autoclavable.

Dimensions: 60 x 64 x 97 mm

Accessories for slide staining box (Order No. 2-1841 ff.)

Additional Information

No.	2-1846
Manufacturer (Brand)	neoLab
EAN	4058072058258
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	gray
Material	Polyoxymethylene (POM)
Wide	60 mm
Height	97 mm
Length	64 mm
Depth	64 mm
Dyeing boxes and racks features	with special chemical resistance
for slide dimensions	76 x 26 mm

