



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® embedding cassettes Macrosette with lid, white, 250 pcs./pack

€124.60
plus VAT &
Shipping

Product Images



Description

Macro embedding cassettes made of high-density polymer for larger sections.
The one-piece lid can be opened several times and closed tightly again.
Specially shaped flow-through openings improve fluid exchange and optimize embedding.
Can be labeled on 3 sides, with large front labeling field.
12 mm high, non-cytotoxic . Solvent resistant.

Color: white

Additional Information

No.	7-0009
Manufacturer (Brand)	neoLab
EAN	4058072026486
Transport temperature	Room temperature
Color	White
Material	Polyoxymethylene (POM)
Wide	40 mm
Height	13 mm
Length	28.5 mm
Depth	28.5 mm

