



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® embedding cassettes for routine embedding, yellow, 500 pcs./pack

€78.70
plus VAT &
Shipping

Product Images



Description

Cassettes made of high-density polymer for embedding or storing samples. Specially shaped flow-through openings improve liquid exchange and optimize embedding. The cassettes can be opened several times and securely resealed.

Solvent resistant.

With large labeling field.

Additional Information

No.	7-0013
Manufacturer (Brand)	neoLab
EAN	4058072026523
Transport temperature	Room temperature
Color	yellow
Material	Polyoxymethylene (POM)
Wide	40 mm
Height	6.5 mm
Length	28.5 mm
Depth	28.5 mm

