



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® staining chamber with transp. lid f. 20 slides

**€513.30**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Suitable not only for routine staining, but also as a humidity chamber for hematology, cytology and microbiology. Base made of black ABS plastic, resistant to many chemicals. Each place is provided with a polymer strip, so that the slides do not slip even when placed at an angle. Bottom drain with stopper for faster emptying. 4 rubber feet for firm stand.

### **Order No. 2-1878**

For 20 slides

**Dimensions:** 380 x 240 x 45 mm

## Additional Information

---

No.	2-1878
Manufacturer (Brand)	neoLab
EAN	4058072058555
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	Black
Material	Acrylonitrile butadiene styrene (ABS)
Wide	240 mm
Height	45 mm
Length	380 mm
Depth	380 mm
Dyeing boxes and racks features	with special chemical resistance with closure/lid
Color closure/ lid	transparent (clear)
Material closure/ lid	Polyethylene terephthalate glycol (PETG)
TBST MAX dyeing chamber	60 °C
TBST MAX closure/ lid	60 °C
TBST MIN Staining chamber	-20 °C
TBST MIN Closure/ lid	-20 °C
for slide dimensions	76 x 26 mm

