



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® staining chamber with black lid f. 20 slides

€513.30
plus VAT &
Shipping

Product Images



Description

Suitable not only for routine staining, but also as a humidity chamber for hematology, cytology and microbiology. Base made of black ABS plastic, resistant to many chemicals. Each place is provided with a polymer strip, so that the slides do not slip even when placed at an angle. Bottom drain with stopper for faster emptying. 4 rubber feet for firm stand.

Order No. 2-1879

For 20 slides

Dimensions: 380 x 240 x 45 mm

Additional Information

No.	2-1879
Manufacturer (Brand)	neoLab
EAN	4058072058562
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	Black
Material	Acrylonitrile butadiene styrene (ABS)
Wide	240 mm
Height	45 mm
Length	380 mm
Depth	380 mm
Dyeing boxes and racks features	with special chemical resistance with closure/lid
Color closure/ lid	black
Capacity	20
Material closure/ lid	Acrylonitrile butadiene styrene (ABS)
TBST MAX dyeing chamber	80 °C
TBST MAX closure/ lid	80 °C
TBST MIN Staining chamber	-80 °C
TBST MIN Closure/ lid	-80 °C
for slide dimensions	76 x 26 mm

