



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® staining box, according to Hellendahl, for 16 slides 76 x 26 mm, soda-lime glass

€272.50
plus VAT &
Shipping

Product Images



Description

Staining trough, Hellendahl type, straight, 90 x 90 x 40 mm Microscopy is a method of viewing objects that are beyond the resolution of the human eye. Objects to be examined are magnified to such an extent that they can be visualized. For sample preparation, the objects to be examined are cut into very thin slices. These slices are then placed on a slide and covered with a coverslip and appropriate solutions. To view the sections under the microscope, the sections on the slide must be stained in a staining trough filled with an appropriate staining solution. Special features and advantages: Can be covered with lid during staining to reduce solvent loss due to spillage or evaporation Staining of up to 10 slides with standard dimensions 76 x 26 mm possible at the same time Application: Staining, Manual staining procedures, Staining of histological and cytological specimens

Additional Information

No.	SD-0978
Reg. trademark	DURAN ®
Manufacturer (Brand)	DURAN
VGKL number	192331400
EAN	4058072346881
Manufacturer EAN	4032051023215
Gross weight	0.452
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.41200 kg
Material	Lime soda glass
Country of origin	Czech Republic
sterile	No
Packaging volume	0.00613300 cc
Customs tariff number	70179000

