



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hecht glass cylinder with knob lid suitable for 2-1971, 120 x 170 mm high

**Price on
request
plus VAT &
Shipping**

Product Images

Description

Glass cylinder with button lid

Suitable for staining slides.

To be used as staining chamber.

Suitable for stainless steel slide racks (Order No. 2-1971/72)

Additional Information

No.	WQ-0877
Manufacturer (Brand)	Hecht
old neoLab article no.	2-1977
EAN	4058072443610
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	120 mm
Color	transparent (clear)
Material	Soda-lime glass (normal glass) 1
Height	170 mm
Dyeing boxes and racks features	with closure/lid
Color closure/ lid	transparent (crystal clear)
Capacity	24
Material closure/ lid	Soda-lime glass (normal glass)
Closure / lid properties	with handle / knob
for slide dimensions	76 x 26 mm

