




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Slide staining box, blue

€24.10
plus VAT &
Shipping

Product Images



Description

EasyDip™, a novel system for staining sections.

The acetal polymer staining box is unbreakable and resistant to most staining reagents (except phenol).

The wide base provides great stability.

With a retaining clip (included), the individual boxes can be loosely fitted together, e.g. to maintain a specific order.

The interior has a volume of only 80 ml.

Temperature resistant from -170°C to +121°C, autoclavable.

Dimensions: 64 x 76 x 92 mm

Accessories: Staining rack (Order No. 2-1846)

Additional Information

No.	2-1841
Manufacturer (Brand)	neoLab
EAN	4058072058203
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue
Material	Polyoxymethylene (POM)
Wide	64 mm
Height	92 mm
Length	76 mm
Depth	76 mm
Dyeing boxes and racks features	stackable/collapsible with special chemical resistance with closure/lid
Color closure/ lid	blue
Material closure/ lid	Polyoxymethylene (POM)
TBST MAX dyeing chamber	121 °C
TBST MAX closure/ lid	121 °C
TBST MIN Staining chamber	170 °C
TBST MIN Closure/ lid	-170 °C
for slide dimensions	76 x 26 mm

