



neoLab® Slide staining box, yellow

Product Images



Description

EasyDipTM, a novel system for staining sections.

The acetal polymer staining box is unbreakable and resistant to most staining reagents (except phenol).

The wide base provides great stability.

With a retaining clip (included), the individual boxes can be loosely fitted together, e.g. to maintain a specific order.

The interior has a volume of only 80 ml.

Temperature resistant from -170°C to +121°C, autoclavable.

Dimensions: 64 x 76 x 92 mm

Accessories: Staining rack (Order No. 2-1846)

Additional Information

No.	2-1845
Manufacturer (Brand)	neoLab
EAN	4058072058241
Color	gelb
Material	Polyoxymethylen (POM)
Wide	64 mm
Height	92 mm
Length	76 mm
Depth	76 mm
Dyeing boxes and racks features	stapel-/zusammensteckbar mit besonderer chemikalischer Beständigkeit mit Verschluss/Deckel
Color closure/ lid	gelb
Material closure/ lid	Polyoxymethylen (POM)
TBST MAX dyeing chamber	121 °C
TBST MAX closure/ lid	121 °C
TBST MIN Staining chamber	170 °C
TBST MIN Closure/ lid	-170 °C
for slide dimensions	76 x 26 mm

