




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® glass staining boxes with lid

**€17.40**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Staining box with lid** made of glass.

Accessories: Insert for 10 slides (Order No. E-1702)

Wire handle (Order No. E-1703)

## Additional Information

---

No.	E-1701
Manufacturer (Brand)	neoLab
EAN	4058072066611
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Soda-lime glass (normal glass) 1
Wide	90 mm
Height	70 mm
Length	108 mm
Depth	108 mm
Dyeing boxes and racks features	with special chemical resistance with closure/lid
Color closure/ lid	transparent (clear)
Material closure/ lid	Soda-lime glass (normal glass) 1
for slide dimensions	76 x 26 mm

