



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Menzel SUPERFROST slides cut, white field, 50 pcs./pack

**€4.57**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Glass slides with white labeling field  
SUPERFROST quality.

The light opal coating can be written on with any pen and is clear and distinct.  
The coating is resistant to common solvents, it prevents the glasses from sticking together and scratching.

Size of microscope slides: 76 x 26 mm, 1 mm thick

Size of labeling field: 15 mm long x 26 mm wide

### Edges cut

**PU:** 50 pcs./pack

## Additional Information

---

No.	MZ-0003
Manufacturer (Brand)	Menzel
old neoLab article no.	1-6266
EAN	4066292131019
Transport temperature	Room temperature
Material	Borosilicate glass
Wide	26 mm
Color matte edge	white
Length	26 mm
Slide properties	with mat edge with coating
Thickness	1 mm
Dimensions	26 x 76 mm
Width mat edge	15 mm
Color slide	bright white
Edges	cut
Type slides	SUPERFROST

