



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Menzel SUPERFROST slides cut, blue field, 50 pcs./pack

€4.84
plus VAT &
Shipping

Product Images



Description

Glass slides with blue labeling field
SUPERFROST quality.

The light opal coating can be written on with any pen and is clear and distinct.
The coating is resistant to common solvents, it prevents the glasses from sticking together and scratching.

Size of microscope slides: 76 x 26 mm, 1 mm thick

Size of labeling field: 15 mm long x 26 mm wide

Edges cut

PU: 50 pcs./pack

Additional Information

No.	MZ-0004
Manufacturer (Brand)	Menzel
old neoLab article no.	1-6267
EAN	4066292131026
Transport temperature	Room temperature
Material	Borosilicate glass
Wide	26 mm
Color matte edge	blue
Length	26 mm
Slide properties	with mat edge with coating
Thickness	1 mm
Dimensions	26 x 76 mm
Width mat edge	15 mm
Color slide	bright white
Edges	cut
Type slides	SUPERFROST

