




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand Dispensette® S, Analog, 0.1-1.0 ml

€343.10
plus VAT &
Shipping

Product Images



Description

Bottle top dispenser with a wide range of applications, DE-M marked and autoclavable (121 °C).

The ejection valve with safety ball minimizes the risk of accidental splashing of the medium. With swivel screw cap, rotatable valve block (360 °) and new valve system without additional sealing rings.

Accuracy:

$\leq \pm 0.5 \%$ (3-0715 = $\leq \pm 1.0 \%$)

Coefficient of variation:

$\leq \pm 0.1 \%$ (3-0715 = $\leq \pm 0.2 \%$)

Adjustment to "Ex". Based on the nominal volume of the device.

Scope of delivery: polypropylene adapter, telescopic suction tube, installation wrench, quality certificate.

Version Analog

- analog volume display
- quickly adjustable dosing volume
- easy adjustment with adjustment key

Additional Information

No.	BN-1320
Manufacturer (Brand)	Brand
old neoLab article no.	3-0715
EAN	4058072206918
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	1 ml
Volume MIN	0.1 ml
Display (reading, setting)	analog
Dispenser properties	with rotatable valve block (360°) with certification with adapter
Functionality	mechanical
Scope of delivery	BRAND Dispensette S® Analog, bottle top adapter 1 x GL45, 1 x GL32, 1 x GL28, 1 x S40, 1 x telescopic aspiration tube, installation wrench, quality certificate, operating instructions
Material valve spring	Platinum Iridium
Model Bottle Top Dispenser	Fire Dispensette S
Correctness at Vol. MAX	0.006
Division	0.02 ml
Type Bottle Top Dispenser	Dispenser (analog)
Type certificate	Declaration of conformity (standardization)
Coefficient of variation at Vol. MAX	0.002
for thread/DM	GL 45
adjusted to	"ex"

