



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Varispenser 2, Bottle Top Dispenser 0.2 – 2 mL

€317.00
plus VAT &
Shipping

Product Images



Description

The **Eppendorf Varispenser® 2 without backflow valve** is the ideal solution for your laboratory tasks with aliquots of aggressive solutions of alkalis, acids, bases or solvents. Thanks to chemically resistant fluoroplastics, the bottle-top dispensers offer maximum chemical resistance. They are as rugged as they are accurate when dispensing virtually any liquid.

The addition of a pressure compensating valve cap with a Luer cone for mounting a microfilter prevents potential contamination from the supply bottle. **The dispensers are fully autoclavable.** The GL 45 thread used on all sizes fits most conventional bottle threads.

Available for different volume ranges. Suitable for cell culture, microbiology, protein biochemistry, molecular biology, chemistry, food analysis, protein analysis, and more.

Additional Information

No.	VB-1629
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-9110
EAN	4058072417222
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	2 ml
Volume MIN	0.2 ml
Display (reading, setting)	analog
Dispenser properties	with certification
Functionality	mechanical
Model Bottle Top Dispenser	Eppendorf Varispenser2
Division	0.1 ml
Material valve spring	Stainless steel
Correctness at Vol. MAX	0.005
Type Bottle Top Dispenser	Dispenser (analog)
Type certificate	Calibration certificate from the manufacturer
Coefficient of variation at Vol. MAX	0.001
for thread/DM	GL 32 GL 38 GL 45 DIN 25 DIN 28 DIN 40
adjusted to	"ex"

