



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Varispenser 2x, Bottle Top Dispenser 0.2 – 2 mL

€375.00
plus VAT &
Shipping

Product Images



Description

The Eppendorf Varispenser® 2x with backflow valve is the ideal solution for your laboratory tasks with aliquots of aggressive solutions, alkalis, acids, bases or solvents. Thanks to chemically resistant fluoroplastics, the bottle-top dispensers offer maximum chemical resistance. They are as rugged as they are accurate when dispensing virtually any liquid. The addition of a pressure equalization valve cap with a Luer cone for mounting a microfilter prevents potential contamination from the supply bottle. The dispensers are fully autoclavable. The GL 45 thread used on all sizes fits most conventional bottle threads. Available for different volume ranges. Suitable for cell culture, microbiology, protein biochemistry, molecular biology, chemistry, food analysis, protein analysis, and more.

Additional Information

No.	VB-1660
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-9116
EAN	4058072417536
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	2 ml
Volume MIN	0.2 ml
Dispenser properties	with certification with back metering valve
Display (reading, setting)	analog
Functionality	mechanical
Material valve spring	Stainless steel
Model Bottle Top Dispenser	Eppendorf Varispenser2x
Correctness at Vol. MAX	0.005
Division	0.1 ml
Type Bottle Top Dispenser	Dispenser (analog)
Type certificate	Calibration certificate from the manufacturer
Coefficient of variation at Vol. MAX	0.001
for thread/DM	GL 32 GL 38 GL 45 DIN 25 DIN 28 DIN 40
adjusted to	"ex"

