



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson DISPENSMAN 10 ML (1-10ml)

Special Price

€289.00 was

€373.00

**plus VAT &
Shipping**

Product Images



Description

DISPENSMAN™ – precise and safe bottle-top dispenser for many common laboratory fluids DISPENSMAN™ features an innovative dispensing tap with three position settings. It allows for controlled dispensing of fluids, lossless flushing of the system with fluid, and removal of air in the system. It also has an anti-drip safety feature that returns fluid to the bottle and does not require a protective cap. Volume changes can be adjusted quickly and easily with one hand by pressing to unlock the regulator and sliding it to the desired volume – keeping the bottle securely closed, minimizing the risk of contamination. DISPENSMAN is easily calibrated by the user and is suitable for most liquids, including many acids, bases and solvents. The system is autoclavable and requires no consumables such as pipette tips. Five DISPENSMAN models cover different volume ranges from 0.25 ml to 50 mL. All meet ISO 8655 specifications, and a variety of adapters ensure the DISPENSMAN can be used with most common bottle types.

Additional Information

No.	CF-0034
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4107
EAN	4058072122706
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	10 ml
Volume MIN	1 ml
Dispenser properties	with rotatable valve block (360°) with reverse dosing valve with adapter
Display (reading, setting)	analog
Functionality	mechanical
Material valve spring	Stainless steel
Model Bottle Top Dispenser	Dispensman 10 ml
Division	0.2 ml
Type Bottle Top Dispenser	Dispenser (analog)
for thread/DM	GL 32 GL 45
adjusted to	"ex"

