



Hirschmann® Ceramus® Dispenser 1.0-5.0 ml

Product Images



Description

High-precision bottle dispenser with a wide range of applications for dispensing almost all liquids.

Resistant special piston made of high-purity ceramic: no abrasion, no swelling and no thermal problems of the piston.

High-quality materials guarantee extremely high chemical resistance.

- user-friendly feedback system for fast venting without reagent loss
- fast and precise volume adjustment via locking adjustment ring
- easy handling during serial dosing due to ergonomic handle and excellent gliding properties of the dosing piston
- Combination of self-locking adjustment and permanently adjusted volume steps, no recalibration necessary
- conformity certified (according to annex 12 of the calibration regulations)
- free rotatability of the device
- guided discharge cannula prevents "snapping effect" of the dosing media
- Valves easy to clean and replace
- sterilizable up to 121°C, also completely assembled

Scope of delivery

Complete with screw thread A 45 mm, adapter for A 32, A 38 and S 40, 1 suction hose with union nut, 1 ejection unit, 1 assembly key.

With individual serial number and calibration certificate.

Additional Information

No.	NH-1312
Manufacturer (Brand)	Hirschmann
old neoLab article no.	3-0602
EAN	4058072336547
Volume MAX	5 ml
Volume MIN	1 ml
Display (reading, setting)	analog
Dispenser properties	mit Zertifizierung mit Rückdosierventil
Functionality	mechanisch
Model Bottle Top Dispenser	Hirschmann-ceramus
Division	0.1 ml
Scope of delivery	Hirschmann-ceramus®-Dispenser, Flaschenaufsatzadapter 1 x GL 32, 1 x S40, 1 x Ansaugschlauch mit Überwurfmutter, 1 x Ausstoßeinheit, 1 x Montageschlüssel, Bedienungsanleitung
Material valve spring	Edelstahl
Correctness at Vol. MAX	0.005
Type Bottle Top Dispenser	Dispenser (analog)
Type certificate	Kalibrierzertifikat vom Hersteller Konformitätserklärung (Normierung)
Coefficient of variation at Vol. MAX	0.002
for thread/DM	GL 32 GL 45 S 40
adjusted to	"ex"

