



## Hirschmann® Ceramus® Dispenser 10.0-60.0 ml

### Product Images

---



## Description

---

High-precision bottle dispenser with a wide range of applications for dispensing almost all liquids.

Resistant special piston made of high-purity ceramic: no abrasion, no swelling and no thermal problems of the piston.

High-quality materials guarantee extremely high chemical resistance.

- user-friendly feedback system for fast venting without reagent loss
- fast and precise volume adjustment via locking adjustment ring
- easy handling during serial dosing due to ergonomic handle and excellent gliding properties of the dosing piston
- Combination of self-locking adjustment and permanently adjusted volume steps, no recalibration necessary
- conformity certified (according to annex 12 of the calibration regulations)
- free rotatability of the device
- guided discharge cannula prevents "snapping effect" of the dosing media
- Valves easy to clean and replace
- sterilizable up to 121°C, also completely assembled

### Scope of delivery

Complete with screw thread A 45 mm, adapter for A 32, A 38 and S 40, 1 suction hose with union nut, 1 ejection unit, 1 assembly key.

With individual serial number and calibration certificate.

## Additional Information

|                                      |                                                                                                                                                                                |
|--------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| No.                                  | NH-1328                                                                                                                                                                        |
| Manufacturer (Brand)                 | Hirschmann                                                                                                                                                                     |
| old neoLab article no.               | 3-0605                                                                                                                                                                         |
| EAN                                  | 4058072336660                                                                                                                                                                  |
| Volume MAX                           | 60 ml                                                                                                                                                                          |
| Volume MIN                           | 10 ml                                                                                                                                                                          |
| Display (reading, setting)           | analog                                                                                                                                                                         |
| Dispenser properties                 | mit Zertifizierung mit Rückdosierventil                                                                                                                                        |
| Functionality                        | mechanisch                                                                                                                                                                     |
| Model Bottle Top Dispenser           | Hirschmann-ceramus                                                                                                                                                             |
| Division                             | 1 ml                                                                                                                                                                           |
| Scope of delivery                    | Hirschmann-ceramus®-Dispenser, Flaschenaufsatzadapter 1 x GL 32, 1 x S40, 1 x Ansaugschlauch mit Überwurfmutter, 1 x Ausstoßeinheit, 1 x Montageschlüssel, Bedienungsanleitung |
| Material valve spring                | Edelstahl                                                                                                                                                                      |
| Correctness at Vol. MAX              | 0.005                                                                                                                                                                          |
| Type Bottle Top Dispenser            | Dispenser (analog)                                                                                                                                                             |
| Type certificate                     | Kalibrierzertifikat vom Hersteller Konformitätserklärung (Normierung)                                                                                                          |
| Coefficient of variation at Vol. MAX | 0.002                                                                                                                                                                          |
| for thread/DM                        | GL 32 GL 45 S 40                                                                                                                                                               |
| adjusted to                          | "ex"                                                                                                                                                                           |

