




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hirschmann® Dispenser fixed volume 1 ml

€347.00
plus VAT &
Shipping

Product Images



Description

High-precision bottle dosing device for reproducible dosing of liquids. Very simple and safe volume setting by turning the scale ring (Order No. 3-0931 to 3-0936).

Steam sterilizable up to 121°C, with high-purity ceramic piston Al₂O₃. The ceramic piston guarantees: no abrasion, no swelling, no thermal problems. Suction and discharge tubes made of PTFE.

Complete with screw thread A 32 and 3 adapters S 40, A 28 and A 45.

With quality certificate, conformity certified.

Accessories: Square bottle, brown (Order No. 3-0941 and 3-0942)

Dispenser bottle, round, brown (Order No. 3-0943)

Additional Information

No.	NH-1243
Manufacturer (Brand)	Hirschmann
old neoLab article no.	3-1500
EAN	4058072336455
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	1 ml
Volume MIN	1 ml
Display (reading, setting)	analog
Dispenser properties	with fixed volume with certification with backflow valve
Functionality	mechanical
Scope of delivery	Hirschmann-ceramus® classic bottle-top dispenser fixed or variable, bottle adapter 1 x GL45, 1 x GL32, 1 x GL28, 1 x S40, quality certificate, operating instructions
Material valve spring	Stainless steel
Model Bottle Top Dispenser	Hirschmann-ceramus classic
Correctness at Vol. MAX	0.006
Division	0.05 ml
Type Bottle Top Dispenser	Dispenser (analog)
Type certificate	Declaration of conformity (standardization) calibration certificate from manufacturer
Coefficient of variation at Vol. MAX	0.002
for thread/DM	GL 32 GL 45 S 40
adjusted to	"ex"

