



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hirschmann® EM dispenser PP 2-10 ml

€268.00
plus VAT &
Shipping

Product Images



Description

Bottle dosing unit with precision-ground glass bulb for problem-free aqueous media as well as for organic media that do not attack polypropylene. With variable volume adjustment in appropriate steps. Sterilizable up to 121°C.

Accuracy $\pm 0.6\%$

Coefficient of variation <math>< 0.2\%</math>

Scope of delivery

Complete with screw thread A 32 mm, adapter for A 28, S 40 and A 45, 1 suction hose with union nut, 1 ejection unit, with individual serial number and weighing record. Certificate of conformity.

Additional Information

No.	NH-1315
Manufacturer (Brand)	Hirschmann
old neoLab article no.	3-1420
EAN	4058072336561
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	10 ml
Volume MIN	2 ml
Display (reading, setting)	analog
Dispenser properties	with certification with back metering valve
Functionality	mechanical
Model Bottle Top Dispenser	neoLab-EM-dispenser
Division	0.25 ml
Scope of delivery	Hirschmann® EM-dispenser PP, bottle top adapter 1 x GL45, 1 x GL28, 1 x S40, 1 x suction hose with union nut, installation wrench, calibration certificate, operating instructions
Material valve spring	Platinum Iridium
Correctness at Vol. MAX	0.006
Type Bottle Top Dispenser	Dispenser (analog)
Type certificate	Declaration of conformity (standardization)
Coefficient of variation at Vol. MAX	0.002
for thread/DM	GL 32 GL 45 S 40
adjusted to	"ex"

