



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hirschmann® Bottle-top dispenser EM- dispenser PP, 10 – 60 ml, graduation 1 ml

€387.00
plus VAT &
Shipping

Product Images



Description

Bottle dosing unit with precision-ground glass bulb for problem-free aqueous media as well as for organic media that do not attack polypropylene. With variable volume adjustment in appropriate steps. Sterilizable up to 121°C.

Accuracy $< \pm 0.6\%$

Coefficient of variation $< 0.2\%$

Scope of delivery

Complete with screw thread A 32 mm, adapter for A 28, S 40 and A 45, 1 suction hose with union nut, 1 ejection unit, with individual serial number and weighing record. Certificate of conformity.

Additional Information

| | |
|--------------------------------------|--|
| No. | NH-1326 |
| Manufacturer (Brand) | Hirschmann |
| old neoLab article no. | 3-1421 |
| EAN | 4058072336646 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Volume MAX | 60 ml |
| Volume MIN | 10 ml |
| Display (reading, setting) | analog |
| Dispenser properties | with certification with back metering valve |
| Functionality | mechanical |
| Scope of delivery | Hirschmann® EM-dispenser PP, bottle top adapter 1 x GL45, 1 x GL28, 1 x S40, 1 x suction hose with union nut, installation wrench, calibration certificate, operating instructions |
| Material valve spring | Platinum Iridium |
| Model Bottle Top Dispenser | neoLab-EM-dispenser |
| Correctness at Vol. MAX | 0.006 |
| Division | 1 ml |
| Type Bottle Top Dispenser | Dispenser (analog) |
| Type certificate | Declaration of conformity (standardization) |
| Coefficient of variation at Vol. MAX | 0.002 |
| for thread/DM | GL 32 GL 45 S 40 |
| adjusted to | "ex" |

