



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Calibrex™ universal 520 adjustable bottle top dispenser 0.25-2 ml

€373.36
plus VAT &
Shipping

Product Images



Description

Bottle Top Dispenser Calibrex™ universal 520

The simple but ingenious design of the bottle top dispenser allows the user to handle it comfortably. The volume is adjusted by rotation. The dispensers fit a variety of reagent containers and bottles, ensuring safe liquid dispensing under state-of-the-art dispensing technologies. The dispensers offer universal chemical compatibility.

Three models, cover volume ranges from 0.25 to 10 ml.

- Chemically resistant materials
- Long-term stability of performance
- Optimal user comfort
- Robust simple construction
- Fast volume adjustment
- Easy recalibration
- 360° rotation possible
- Fully assembled autoclavable at 121°C
- QR code for chemical resistance table
- Two-year warranty.

Optional (part no. 7-7453 to 7-7455)

Fixed models with a user-defined volume are also available. The volume selected depends on the respective volume ranges.

Additional Information

No.	SR-0178
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7450
EAN	4058072366964
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	2 ml
Volume MIN	0.25 ml
Display (reading, setting)	analog
Dispenser properties	with fixed volume
Functionality	mechanical
Model Bottle Top Dispenser	Socorex-Calibrex universal TM 520
Scope of delivery	Dispenser with 32 mm standard thread, 3 adapters (PP) : 28, 40, 45 mm, extension tube (300 mm)
Correctness at Vol. MAX	0.005
Type Bottle Top Dispenser	Dispenser (analog)
Coefficient of variation at Vol. MAX	0.001
for thread/DM	GL 25 GL 28 GL 32

