



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex Calibrex™ universal 520 adjustable bottle top dispenser 1-10 ml

**€373.36**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### **Bottle Top Dispenser Calibrex™ universal 520**

The simple but ingenious design of the bottle top dispenser allows the user to handle it comfortably. The volume is adjusted by rotation. The dispensers fit a variety of reagent containers and bottles, ensuring safe liquid dispensing under state-of-the-art dispensing technologies. The dispensers offer universal chemical compatibility.

Three models, cover volume ranges from 0.25 to 10 ml.

- Chemically resistant materials
- Long-term stability of performance
- Optimal user comfort
- Robust simple construction
- Fast volume adjustment
- Easy recalibration
- 360° rotation possible
- Fully assembled autoclavable at 121°C
- QR code for chemical resistance table
- Two-year warranty.

Optional (part no. 7-7453 to 7-7455)

Fixed models with a user-defined volume are also available. The volume selected depends on the respective volume ranges.

## Additional Information

No.	SR-0180
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7452
EAN	4058072366988
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	10 ml
Volume MIN	1 ml
Display (reading, setting)	analog
Dispenser properties	with fixed volume
Functionality	mechanical
Model Bottle Top Dispenser	Socorex-Calibrex universal <sup>TM</sup> 520
Scope of delivery	Dispenser with 32 mm standard thread, 3 adapters (PP) : 28, 40, 45 mm, extension tube (300 mm)
Correctness at Vol. MAX	0.005
Type Bottle Top Dispenser	Dispenser (analog)
Coefficient of variation at Vol. MAX	0.001
for thread/DM	GL 25 GL 28 GL 32

