



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Acurex™ compact 501, Compact Dispenser 0.2-2 ml/Reservoir 2000 ml

€366.08
plus VAT &
Shipping

Product Images



Description

Acurex™ compact compact dispenser for safe precise dispensing of reagents.

Thanks to its small space requirement, the ^{Acurex™} 501 is especially suitable for setting in refrigerators. Can also be used in a water bath.

Volume range from 0.2 to 30 ml.

- Fully protected dispensing mechanism
- Retractable graduated guide rod for reduced space requirements
- Amber glass for optimal light protection
- Best chemical resistance of all parts in contact with liquid, metal-free
- Choice of four reservoir sizes
- Fully assembled autoclavable at 121° C

Two-year warranty.

Additional Information

No.	SR-0171
Manufacturer (Brand)	Socorex
old neolab article no.	7-7442
EAN	4058072366896
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	2 ml
Volume MIN	0.2 ml
Display (reading, setting)	analog
Dispenser properties	with storage vessel
Functionality	mechanical
Model Bottle Top Dispenser	Socorex-Acurex™501
Division	0.1 ml
Volume storage vessel	2000.000000 ml
Correctness at Vol. MAX	0.005
Type Bottle Top Dispenser	Dispenser (analog)
Coefficient of variation at Vol. MAX	0.001

