



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex Calibrex™ organo 525 Bottle Top Dispenser, with Glass Flask, 1-10 ml

**€384.28**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### **Bottle-top dispenser Calibrex™ organo 525**

This superior instrument line guarantees safe as well as flexible distribution of fluids, fitting a variety of reagent bottles. Excellent chemical resistance, robustness, easy handling and maintenance stand for top performance in the laboratory.

Optionally available with fluid control system with stopcock (part no. 7-7487 to 7-7493), which allows easy venting without loss of reagents.

Seven models cover a volume range from 0.1 to 100 ml.

The Calibrex™ organo has a ground glass bulb and is particularly suitable for organic and non-crystallizing solvents.

For detailed information on chemical resistance, please refer to the chemical resistance table.

- Fluid control system with stopcock (optional)
- Permanent visibility of fluid flow
- Integrated calibration system
- Freely rotatable 360° on the bottle
- Volume adjustment by means of slide regulator (optionally replaceable; screw regulator)
- Long-term stability of performance
- QR code imprint for access to chemical resistance chart.
- Fully assembled autoclavable at 121°C.
- Two-year warranty.

## Additional Information

|                                      |  |
|--------------------------------------|--|
| No.                                  | SR-0207  |
| Manufacturer (Brand)                 | Socorex  |
| old neolab article no.               | 7-7483   |
| EAN                                  | 4058072367251  |
| Transport temperature                | Room temperature   |
| autoclavable at 121°C                | Yes  |
| Volume MAX                           | 10 ml  |
| Volume MIN                           | 1 ml   |
| Display (reading, setting)           | analog   |
| Dispenser properties                 | with rotatable valve block (360°) with adapter   |
| Functionality                        | mechanical   |
| Model Bottle Top Dispenser           | Socorex-Calibrex™525   |
| Scope of delivery                    | Dispenser, feeding tube, dosing cannula, manufacturer's quality certificate, operating instructions. Standard thread 45 mm, adapters in different sizes. |
| Correctness at Vol. MAX              | 0.006  |
| Type Bottle Top Dispenser            | Dispenser (analog)   |
| Coefficient of variation at Vol. MAX | 0.001  |
| for thread/DM                        | GL 28 GL 32 GL 40  |

