



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex Calibrex™ organo 525 Bottle Top Dispenser, with glass bulb, 2.5-25 ml

**€491.40**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### **Bottle-top dispenser Calibrex™ organo 525**

This superior instrument line guarantees safe as well as flexible distribution of fluids, fitting a variety of reagent bottles. Excellent chemical resistance, robustness, easy handling and maintenance stand for top performance in the laboratory.

Optionally available with fluid control system with stopcock (part no. 7-7487 to 7-7493), which allows easy venting without loss of reagents.

Seven models cover a volume range from 0.1 to 100 ml.

The Calibrex™ organo has a ground glass bulb and is particularly suitable for organic and non-crystallizing solvents.

For detailed information on chemical resistance, please refer to the chemical resistance table.

- Fluid control system with stopcock (optional)
- Permanent visibility of fluid flow
- Integrated calibration system
- Freely rotatable 360° on the bottle
- Volume adjustment by means of slide regulator (optionally replaceable; screw regulator)
- Long-term stability of performance
- QR code imprint for access to chemical resistance chart.
- Fully assembled autoclavable at 121°C.
- Two-year warranty.

## Additional Information

No.	SR-0208
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7484
EAN	4058072367268
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	25 ml
Volume MIN	2.5 ml
Display (reading, setting)	analog
Dispenser properties	with rotatable valve block (360°) with adapter
Functionality	mechanical
Scope of delivery	Dispenser, feeding tube, dosing cannula, manufacturer's quality certificate, operating instructions. Standard thread 45 mm, adapters in different sizes.
Model Bottle Top Dispenser	Socorex-Calibrex™525
Correctness at Vol. MAX	0.006
Type Bottle Top Dispenser	Dispenser (analog)
Coefficient of variation at Vol. MAX	0.001
for thread/DM	GL 28 GL 32 GL 40

