



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Calibrex™ organo 525 bottle top dispenser, with glass bulb, stopcock, 2.5-25 ml

€548.60
plus VAT &
Shipping

Product Images



Description

Bottle-top dispenser Calibrex™ organo 525

This superior instrument line guarantees safe as well as flexible distribution of fluids, fitting a variety of reagent bottles. Excellent chemical resistance, robustness, easy handling and maintenance stand for top performance in the laboratory.

Optionally available with fluid control system with stopcock (part no. 7-7487 to 7-7493), which allows easy venting without loss of reagents.

Seven models cover a volume range from 0.1 to 100 ml.

The Calibrex™ organo has a ground glass bulb and is particularly suitable for organic and non-crystallizing solvents.

For detailed information on chemical resistance, please refer to the chemical resistance table.

- Fluid control system with stopcock (optional)
- Permanent visibility of fluid flow
- Integrated calibration system
- Freely rotatable 360° on the bottle
- Volume adjustment by means of slide regulator (optionally replaceable; screw regulator)
- Long-term stability of performance
- QR code imprint for access to chemical resistance chart.
- Fully assembled autoclavable at 121°C.
- Two-year warranty.

Additional Information

| | |
|--------------------------------------|--|
| No. | SR-0215 |
| Manufacturer (Brand) | Socorex |
| old neoLab article no. | 7-7491 |
| EAN | 4058072367336 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Volume MAX | 25 ml |
| Volume MIN | 2.5 ml |
| Display (reading, setting) | analog |
| Dispenser properties | with rotatable valve block (360°) with adapter |
| Functionality | mechanical |
| Scope of delivery | Dispenser, feeding tube, dosing cannula, manufacturer's quality certificate, operating instructions. Standard thread 45 mm, adapters in different sizes. |
| Model Bottle Top Dispenser | Socorex-Calibrex™525 |
| Correctness at Vol. MAX | 0.006 |
| Type Bottle Top Dispenser | Dispenser (analog) |
| Coefficient of variation at Vol. MAX | 0.001 |
| for thread/DM | GL 28 GL 32 GL 40 |

