



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## **sunlab® Bottle Top Dispenser, mechanical, 2.5-25 ml**

**€240.00**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

The Sunlab® bottle-top dispenser is suitable for a wide range of dispensing tasks in the laboratory. Due to its high chemical resistance, most acids, alkalis and solvents can be easily dispensed. All parts that come into contact with liquids are made of resistant plastics (BSG, PTFE, FEP). The dispenser is easy to clean and fully autoclavable. Various adapters for different thread types and sizes included in the scope of delivery make the Sunlab® Bottle Top Dispenser suitable for a wide range of bottles and vessels.

## Additional Information

No.	D-8625
Manufacturer (Brand)	sunlab
EAN	4058072086404
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	25 ml
Volume MIN	2.5 ml
Display (reading, setting)	analog
Spare part	D-8651 D-8652
Functionality	mechanical
Scope of delivery	Sunlab® Bottle Top Dispenser, Bottle Top Adapter 1 x S40, 1 x GL32, 1 x GL28, 1 x GL25, 1 x Ejection Cannula, 1 x Aspiration Tube, Operating Instructions
Material valve spring	Platinum/Iridium
Model Bottle Top Dispenser	Sunlab Bottle Top Dispenser
Correctness at Vol. MAX	0.005
Sunlab® Model number	Bottle top dispenser SUI605
Division	2.5 ml
Type Bottle Top Dispenser	Dispenser (analog)
Coefficient of variation at Vol. MAX	0.001
for thread/DM	GL 45 GL 32 GL 28 GL 25 S 40
adjusted to	"ex"

