



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Vitlab neoSpenser, 0.5 - 2.5 ml, subdivision 0.05 ml

Special Price

€125.00 was

€205.00

plus VAT &
Shipping

Product Images



Description

neoSpenser, the low-cost bottle-top dispenser.

- easy and precise volume adjustment
- resistant to most acids, alkalis and solvents
- direct displacement piston with PFA seal reduces crystallization of liquids
- can be disassembled without tools
- high-quality materials: PTFE, platinum-iridium, borosilicate glass
- fully autoclavable at 121°C, 2 bar, max. 20 min.
- conformity certified according to DIN12600
- individual quality certificate

Systematic measurement deviation R %:

$\leq + 0.6$ or 0.5

Random measurement deviation V %: ≤ 0.1

Scope of delivery

Adapter for GL 28, 45 and S 40 with nominal volume 2.5/5/10 ml

Adapter for GL 32, 38 and S 40 with nominal volume 25 ml

Complete with ejection cannula, aspiration tube, 3 adapters and tool.

Accessories: Square bottle, brown (Order No. 3-0941 and 3-0942)

Dispenser bottle, round, brown (Order No. 3-0943)

Additional Information

No.	LV-1820
Manufacturer (Brand)	Vitlab
old neoLab article no.	7-0086
EAN	4066292052963
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	2.5 mL
Volume MIN	0.5 ml
Display (reading, setting)	analog
Dispenser properties	with certification
Functionality	mechanical
Scope of delivery	Bottle top dispenser, bottle adapter 1 x GL45, 1 x GL32, 1 x GL28, 1 x S40, operating instructions
Material valve spring	Platinum Iridium
Model Bottle Top Dispenser	neoSpenser
Correctness at Vol. MAX	0.006
Division	0.05 ml
Type Bottle Top Dispenser	Dispenser (analog)
Type certificate	Declaration of conformity (standardization) calibration certificate from manufacturer
Coefficient of variation at Vol. MAX	0.002
for thread/DM	GL 32 GL 45 S 40
adjusted to	"ex"

