



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Brand HandyStep® S, incl. shelf holder

**€250.35**  
**plus VAT &**  
**Shipping**

## Product Images



## Description

---

The **HandyStep® S manual multiple dispenser makes serial pipetting even easier**. With just one filling, volumes between 2 and 5 ml can be dispensed with the PD tips in up to 49 pipetting steps without refilling.

- stroke adjustment wheel with 9 positions
- 59 different partial volumes can be set with 10 PD-Tip sizes (incl. 1 and 10 ml tips)
- 2-sided volume table on the back of the device for quick volume determination
- easy tip coupling - PD-Tips are inserted straight from below
- safe tip ejection with only one button - without touching the dispenser tips
- clear display, no battery necessary
- increased chemical resistance due to innovative plastics
- maintenance free, factory adjusted
- best possible results are achieved with PD-Tips from BRAND. Can also be used with Combitips® plus, Combitips®, Repet-Tips, Encode™-Tips and other compatible tips

Ideal for routine and research in diagnostics, molecular biology, environmental analysis, etc. HandyStep® S and PD-Tips work according to the direct displacement principle. This allows liquid media with high viscosity, high density or high vapor pressure to be dispensed with maximum precision.

## Additional Information

No.	BN-1770
Manufacturer (Brand)	Brand
old neoLab article no.	3-0825
EAN	4058072210328
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	light gray
Volume MAX	50000 µl
Volume MIN	0.002 µl
sterile	No
Operation via	Knobs/Buttons
Dispenser properties	with certification with dispensing function ergonomically shaped
Number of channels	1
Energy supply	without power supply
Functionality	mechanical
Model Microliter pipette/ Dispenser	HandyStep S
Division	2 µl
Type Microliter pipette (design)	Direct displacement pipette
Manufacturer	Fire
Scope of delivery	HandyStep® S, shelf holder, 1 x PD-Tip 0.1ml, 1 x PD Tip 1 ml, 1 x PD-Tip 10ml, quality certificate, operating instructions
Type pipette	Microliter pipettes

