



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand multiple dispenser HandyStep touch S, incl. universal power supply and retaining clip

**€1,017.45
plus VAT &
Shipping**

Product Images



Description

HandyStep® touch S

The multiple dispenser from BRAND

The new multiple dispenser from BRAND is the first choice when it comes to efficient dispensing of long series.

HandyStep® touch S with the additional functions Sequential Dispensing, Multi-Aspiration and Titration

At a glance

- Touch operation with intuitive menu guidance and integrated help
- Available as HandyStep® touch or HandyStep® touch S model with extended range of functions
- Languages: English, German, Spanish, French and Chinese
- Variable volume setting from 1.0 µl to 50 ml
- Inductively rechargeable

Enduring

Many hours of continuous operation with just one battery charge

Robust

Long service life due to high-quality materials and easy cleaning/disinfection due to smooth surfaces.

Safe

Automatic size recognition of BRAND PD-Tips II

Efficient

In addition to PD-Tips from BRAND, compatible dispenser tips from other manufacturers can also be used.

Uncomplicated

Automatic tip ejection at the touch of a button avoids contamination.

Convenient

Inductive charging with optional charging stand

Step button

Ergonomically positioned step button for effortless volume pick-up and dispensing

Additional Information

No.	BN-2615
Manufacturer (Brand)	Brand
EAN	4058072217464
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	Black
Material	Polycarbonate (PC) Polypropylene (PP)
Volume MAX	50 ml
Volume MIN	1 µl
sterile	No
Display (reading, setting)	graphic display
Operation via	Touchscreen
Dispenser properties	including batteries/battery with charging station with dispensing function with titration function with certification ergonomically shaped
Number of channels	1
Energy supply	Battery
Functionality	electronic
Number of parts in set	5
Type Microliter pipette (design)	Direct displacement pipette
Spare part	BN-2617 BN-2620 BN-2619
Manufacturer	Fire
Scope of delivery	HandyStep® touch S, DE-M marked, with quality certificate, user manual, quick start guide, wall/shelf mount, Li-battery, universal power supply and USB type C cable. One PD-Tip each size 0.5 ml, 1.25 ml, 2.5 ml, 5 ml and 12.5 ml.
Correctness at Vol. MAX	0.03
Type certificate	Declaration of conformity (standardization)

