




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Eppendorf Dispenser Multipette M4, , for Combitips advanced

**€420.00**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

With variable volume selection, you can choose from up to 100 volumes.

Combitips advanced also allow dispensing in up to 100 increments with just one filling, depending on the volume.

The Multipette®M4 covers volume sizes from 0.1 ml to 50 ml with 9 Combitips advanced.

A volume spectrum from 1 µl to 10 ml.

- simplified volume setting with volume selection wheel
- Volume display with Combitip recognition via display
- Direct displacement principle
- easy Combitip pick-up
- ergonomically shaped operating and filling lever
- Motion sensor switches off the Multipette during breaks in use

## Additional Information

No.	VB-1690
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-2605
EAN	4058072417796
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Volume MAX	50000 µl
Volume MIN	1 µl
Operation via	Knobs/Buttons
Dispenser properties	with certification with dispensing function including batteries/battery ergonomically shaped
Number of channels	1
Energy supply	Battery
Functionality	mechanical
Manufacturer	Eppendorf
Measurement deviation systematic in %	0.3 %
Measurement deviation systematic in µl	30 µl
Measurement deviation random in %	0.3 %
Measurement deviation random in µl	30 µl
Model Microliter pipette/ Dispenser	Multipette M4
Division	1 µl
Type Microliter pipette (design)	Direct displacement pipette
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization)

