



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Brand Package HandyStep touch + Charging Stand Universal Power Supply/Universal Holder, White

**€894.60**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Multiple dispensers are the first choice when it comes to efficient dispensing of long series. For pipetting viscous media, direct displacement pipettes are often the right solution. The Stepper HandyStep® touch from BRAND combines both features and can do even more. With the color touch display, you can conveniently flick back and forth between functions and always have all the important information about your work steps in view. Operation is child's play even when wearing protective gloves. Speaking of easy:

the days of changing tips manually are over. With the HandyStep® touch, all it takes is a light tap on the display and the tip is ejected. The HandyStep® touch multiple dispenser automatically detects the volume of the original BRAND PD-Tips II. This saves time and reduces the risk of errors. For the greatest possible flexibility, however, nothing stands in the way of using tips from many other manufacturers. The HandyStep® touch is also simple and efficient when it comes to charging:

once positioned in the charging stand, the HandyStep® touch battery charges inductively. At a glance - Touch operation with intuitive menu navigation - Languages:

German, English, Spanish, French and Chinese - HandyStep® touch with the functions multi-dispensing, auto-dispensing, pipetting, HandyStep® touch S additionally with the functions sequential dispensing, multi-aspiration and titration - Automatic size recognition of BRAND PD-Tips

II - Open for tips from many other manufacturers, e. g.E.g. Encode™ tips, Repet tips, Combitips®, Combitips® plus, Combitips® advanced (Represented trademarks are property of their respective owners) - Automatic tip ejection - Variable volume setting from 1.0 µl to 50 ml - Easy to clean thanks to smooth surface without edges and grooves - Inductively rechargeable

## Additional Information

No.	BN-2703
Article - USP	Simple: touch operation like on a smartphone Uncomplicated: automatic tip ejection Efficient: open for tips from other manufacturers
Manufacturer (Brand)	Brand
VGKL number	104975201
EAN	4058072218072
Gross weight	1.833
HSNumber	84798997
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	1.30000 kg
Color	White
UN number	3481
Country of origin	Germany
Packing width	0.17000 m
Packing height	0.12000 m
Packing depth	0.25000 m
Packaging volume	5100.000000 ccm
Customs tariff number	84798997
Youtube URL	<a href="https://youtu.be/lISsSaOoM-c">https://youtu.be/lISsSaOoM-c</a>

