




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex Stepper™ 411 repeating pipette, 10–5000 $\mu$ l

**Price on  
request  
plus VAT &  
Shipping**

## Product Images



## Description

With this highly convenient **Stepper™ 411 repetitive pipette** (hand dispenser), you can dispense precisely and obtain flawless results. The unique four-finger pressure function relieves the thumb. Carefully selected materials provide exceptional shock resistance. The compact and simple design of the housing reduces the risk of corrosion and internal contamination. Up to 73 dispensings are possible within a choice of 53 different volumes.

Volume range from 10 to 5000 µl.

- Four-finger operation instead of thumb
- Adjustment knob with clear indication of volume and number of desired doses
- Three color-coded direct displacement syringes cover the complete volume range
- Safety lock at the end of the dispensing process

Two years warranty

## Additional Information

|                                  |                                                                          |
|----------------------------------|--------------------------------------------------------------------------|
| No.                              | SR-0147                                                                  |
| Manufacturer (Brand)             | Socorex                                                                  |
| old neolab article no.           | 7-7414                                                                   |
| EAN                              | 4066292207950                                                            |
| Transport temperature            | Room temperature                                                         |
| autoclavable at 121°C            | Yes                                                                      |
| Color                            | White                                                                    |
| sterile                          | No                                                                       |
| Operation via                    | Knobs/Buttons                                                            |
| Dispenser properties             | with certification with dispensing function ergonomically shaped         |
| Number of channels               | 1                                                                        |
| Energy supply                    | without power supply                                                     |
| Functionality                    | mechanical                                                               |
| Manufacturer                     | Socorex                                                                  |
| Scope of delivery                | Stepper™ 411, 3 color-coded adjustment knobs, adapter, Ecostep™ syringes |
| Type Microliter pipette (design) | Direct displacement pipette                                              |
| Type certificate                 | Sterilization certificate from the manufacturer                          |

