



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Acura® capillary 846 direct displacement micropipette, volume 5-25 µl

€318.24
plus VAT &
Shipping

Product Images



Description

Socorex Acura® capillary 846 **direct displacement micropipettes** are highly efficient instruments with replaceable glass capillaries and ETFE-coated plunger. The ergonomic shape as well as gentle activation of the functions ensure optimal hand comfort. The pipette is used for pipetting dense, viscous, foaming or volatile liquids.

- Five models in the volume range from 1 to 200 µl
- Slim shape, comfortable finger support
- Smooth piston stroke
- Reliable stepwise volume adjustment
- Matching color coding on pipettes, pistons and capillaries
- Stainless steel flask, ETFE-coated

Additional Information

No.	SR-0131
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7398
EAN	4058072366612
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Volume MAX	25 µl
Volume MIN	5 µl
Operation via	Knobs/Buttons
Microliter pipettes features	with color code
Number of channels	1
Functionality	mechanical
Type Microliter pipette (design)	Direct displacement pipette
Manufacturer	Socorex
Measurement deviation systematic in %	0.8 %
Measurement deviation systematic in µl	0.2 µl
Measurement deviation random in %	0.6 %
Measurement deviation random in µl	0.15 µl

