



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® BagPipet® adjustable 0.1 ml, 0.9 ml, 1.0 ml

**€182.85**  
plus VAT &  
Shipping

### Product Images



## Description

---

The ideal complement to the <sup>BagMixer®</sup> !

BagPipet®, an automatic pipette, adjustable to 0.1, 0.9 or 1.0 ml.

## Additional Information

---

No.	2-1395
Manufacturer (Brand)	neoLab
EAN	4058072056476
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	1000 µl
Volume MIN	100 µl
Operation via	Knobs/Buttons
Microliter pipettes features	ergonomically shaped
Number of channels	1
Energy supply	without power supply
Functionality	mechanical
Model Microliter pipette/ Dispenser	BagPipet
Type Microliter pipette (design)	Piston stroke pipette/air cushion pipette
Manufacturer	neoLab
Type pipette	Microliter pipettes

