



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Reference® 2, 0.1 – 2.5 µL, dark gray

€376.00
plus VAT &
Shipping

Product Images



Description

The Reference®2 microliter pipette offers exceptional precision.

With its ergonomic design, variable design for volumes and long service life, it is a reliable partner for demanding work in diagnostics, research and industry.

- Simple and visual volume adjustment with volume fixation
- Volume range from 0.1 µl to 10 ml
- Fully autoclavable
- Single button operation
- Reproducible tip fit on the pipette cone

Each variable Reference up to 1000 µl contains an epT.I.P.S. box with tips of the corresponding size.

Pipettes with a larger final volume (2.5/5.0/10 ml) are supplied with an epT.I.P.S.® sample bag.

Additional Information

No.	VB-1489
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-0430
EAN	4058072416010
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Volume MAX	2.5 µL
Volume MIN	0.1 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with GLP functions with color code ergonomically shaped
Number of channels	1
Energy supply	without power supply
Color code	dark gray
Functionality	mechanical
Model Microliter pipette/ Dispenser	Reference 2 0.1 - 2.5µl
Division	0.002 µl
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Manufacturer	Eppendorf
Measurement deviation systematic in %	1.4 %
Measurement deviation systematic in µl	0.04 µl
Measurement deviation random in %	0.7 %
Measurement deviation random in µl	0.018 µl
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization)

