



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Reference® 2, 10 – 100 µL, yellow

**€376.00**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

The Reference®2 microliter pipette offers exceptional precision.

With its ergonomic design, variable design for volumes and long service life, it is a reliable partner for demanding work in diagnostics, research and industry.

- Simple and visual volume adjustment with volume fixation
- Volume range from 0.1 µl to 10 ml
- Fully autoclavable
- Single button operation
- Reproducible tip fit on the pipette cone

Each variable Reference up to 1000 µl contains an epT.I.P.S. box with tips of the corresponding size.

Pipettes with a larger final volume (2.5/5.0/10 ml) are supplied with an epT.I.P.S.® sample bag.

## Additional Information

No.	VB-1493
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-0434
EAN	4058072187361
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Volume MAX	100 µl
Volume MIN	10 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with GLP functions with color code ergonomically shaped
Number of channels	1
Energy supply	without power supply
Color code	yellow
Functionality	mechanical
Manufacturer	Eppendorf
Measurement deviation systematic in %	0.8 %
Measurement deviation systematic in µl	0.8 µl
Measurement deviation random in %	0.2 %
Measurement deviation random in µl	0.2 µl
Model Microliter pipette/ Dispenser	Reference 2 10 - 100µl
Division	0.1 µl
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Type pipette	Microliter pipettes
Type certificate	Declaration of conformity (standardization)

