



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Research® plus, fix, 20 µL, light gray

€190.00
plus VAT &
Shipping

Product Images



Description

The ultra-light Eppendorf Research plus is one of the most frequently used pipettes in the world. In addition to safety and robustness, its development focused on ergonomic operation.

For fixed volume

- **Applications**
- Forward pipetting
- Reverse pipetting
- Removal of supernatants
- Sample mixing
- Phase extraction
- Volume determination
- Filling of plates, gels and reaction vessels

Additional Information

No.	VB-0952
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-0872
EAN	4058072171384
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Volume MAX	20 µl
Volume MIN	20 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with color code ergonomically shaped
Number of channels	1
Energy supply	without power supply
Color code	light gray
Functionality	mechanical
Division	0.02 µl
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Manufacturer	Eppendorf
Measurement deviation systematic in %	0.8 %
Measurement deviation systematic in µl	0.16 µl
Measurement deviation random in %	0.3 %
Measurement deviation random in µl	0.06 µl
Type pipette	Microliter pipettes
Type certificate	Calibration certificate from manufacturer Declaration of conformity (standardization)

