



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Research® plus, fix, 50 µL, yellow

€190.00
plus VAT &
Shipping

Product Images



Description

The ultra-light Eppendorf Research plus is one of the most frequently used pipettes in the world. In addition to safety and robustness, its development focused on ergonomic operation.

For fixed volume

- **Applications**
- Forward pipetting
- Reverse pipetting
- Removal of supernatants
- Sample mixing
- Phase extraction
- Volume determination
- Filling of plates, gels and reaction vessels

Additional Information

No.	VB-0955
Manufacturer (Brand)	Eppendorf
old neolab article no.	E-0875
EAN	4058072411060
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Volume MAX	50 µl
Volume MIN	50 µl
sterile	No
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with color code ergonomically shaped
Number of channels	1
Energy supply	without power supply
Color code	yellow
Functionality	mechanical
Manufacturer	Eppendorf
Measurement deviation systematic in %	0.7 %
Measurement deviation systematic in µl	0.35 µl
Measurement deviation random in %	0.3 %
Measurement deviation random in µl	0.15 µl
Division	0.1 µl
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Type pipette	Microliter pipettes
Type certificate	Calibration certificate from manufacturer Declaration of conformity (standardization)

