



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Xplorer single-channel pipette, , 15-300 µl

**€700.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Single-channel pipette, type Variable

The new electronic pipette Eppendorf<sup>Xplorer®</sup> sets new standards in terms of simplicity, precision and reproducibility. It has been specially developed for high professional requirements.

This is largely achieved by the new intuitive operating concept and the ergonomic design based on the successful Eppendorf PhysioCare<sup>Concept®</sup>. These features set new standards in terms of simplicity, precision and reproducibility. No more wasted time due to complicated programming or rigid procedures. Instead, precisely adjustable parameters, maximally reproducible results, fatigue-free work and always full control over the pipette's processes.

The large and ergonomically positioned multifunction rocker with its "up is up and down is down"<sup>™</sup> principle guarantees uniform and reproducible aspiration and dispensing of liquids, and at the same time serves for color coding. Volume and speed can be set precisely and reproducibly.

- Dial: all functions visible and selectable at a glance
- multilingual menu navigation
- tidy color display
- every adjustable parameter is visible on the display, no submenus
- innovative ejector: pipette automatically moves to zero position after tip ejection
- Quick release: lower part easy to remove and maintain
- spring-loaded tip cone: optimum tightness with minimum attachment force
- Volume & speed can be adjusted precisely and reproducibly
- Adaptation to individual environments, for example, by adjusting for a particular liquid density or a different geographic altitude
- Channel display: Use the pipette in the same orientation - all the time
- remove individual channels: Gain flexibility and safety
- intuitive operating concept

## Additional Information

No.	VB-1109
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-1842
EAN	4058072412296
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	light gray
Volume MAX	300 µl
Volume MIN	15 µl
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with mains/charging adapter with color code with dispensing function including batteries/battery ergonomically shaped IVD compliant
Number of channels	1
Energy supply	Battery Mains
Color code	orange
Functionality	electronic
Model Microliter pipette/ Dispenser	Xplorer 300µl
Division	0.3 µl
Type Microliter pipette (design)	Piston stroke pipette air cushion pipette
Manufacturer	Eppendorf
Measurement deviation systematic in %	0.6 %
Measurement deviation systematic in µl	1.8 µl
Measurement deviation random in %	0.2 %
Measurement deviation random in µl	0.6 µl
Type pipette	Microliter pipettes
Type certificate	Declaration of Conformity (standardization) IVD Directive 98/79/EC

