



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson PIPETMAN® G P2G microliter pipette, metal ejector, volume 0.2-2 µL

€293.00
plus VAT &
Shipping

Product Images



Description

Single-channel pipette, type Variable

Pipetman® G Premium comfort in legendary Gilson quality

- lowest effort during pipetting
- Reduction of tip ejection forces
- Tip ejector with clip mechanism
- Tip ejector also available in plastic

All volume ranges available:

8 PIPETMAN G models comfortably cover the range from 0.2 µl to 10 ml.

Accessories: Carousel stand for 7 pipettes (Order No. 7-4307)

Additional Information

No.	CF-0169
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4415
EAN	4058072220426
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Volume MAX	2 µl
Volume MIN	0.2 µl
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with GLP functions with color code ergonomically shaped IVD compliant
Number of channels	1
Energy supply	without power supply
Color code	orange
Functionality	mechanical
Model Microliter pipette/ Dispenser	Pipetman G P2G
Division	0.002 µl
Designation tip (Gilson)	D10 DL10 DF10ST DFL10ST
Type Microliter pipette (design)	Air cushion pipette
Manufacturer	Gilson
Measurement deviation systematic in %	1.5 %
Measurement deviation systematic in µl	0.03 µl
Measurement deviation random in %	2 %
Measurement deviation random in µl	0.04 µl
Type pipette	Microliter pipettes
Type certificate	Declaration of Conformity (standardization) IVD Directive 98/79/EC

