



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Gilson Pipetman® P100G, 10–100 µl, plastic ejector

**€222.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Pipetman® G with plastic ejector Premium comfort in legendary Gilson quality.

Minimal force required for pipetting.

Reduction of tip ejection forces.

Tip ejector with clip mechanism.

## Additional Information

No.	CF-0880
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4513
EAN	4066292499249
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Volume MAX	10 µl
Volume MIN	100 µl
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with GLP functions with color code ergonomically shaped IVD compliant
Number of channels	1
Designation tip (Gilson)	D200 DF200ST
Energy supply	without power supply
Color code	pink
Functionality	mechanical
Manufacturer	Gilson
Measurement deviation systematic in %	1 %
Measurement deviation systematic in µl	0.8 µl
Measurement deviation random in %	0.3 %
Measurement deviation random in µl	0.3 µl
Model Microliter pipette/ Dispenser	Pipetman P100G with plastic ejector
Division	0.1 µl
Type Microliter pipette (design)	Air cushion pipette
Type pipette	Microliter pipettes
Type certificate	Declaration of Conformity (standardization) IVD Directive 98/79/EC

