



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## **Gilson PIPETMAN® L P20L microliter pipette, metal ejector, volume 2-20 µL**

**€327.00**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

Single-channel pipette, variable type

More productivity more convenience

- ergonomic design
- light in weight
- with metal ejector
- 3 tip ejector positions - ideal for left- and right-handed users

GLP Features

- Volume lock/unlock mechanism
- unique 2D identification code
- traceable serial number
- color coded push buttons

## Additional Information

No.	CF-0616
Manufacturer (Brand)	Gilson
old neoLab article no.	7-3402
EAN	4058072222888
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Volume MAX	20 µl
Volume MIN	2 µl
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with GLP functions with color code ergonomically shaped IVD compliant
Number of channels	1
Energy supply	without power supply
Color code	yellow
Functionality	mechanical
Model Microliter pipette/ Dispenser	Pipetman L P20L
Division	0.02 µl
Designation tip (Gilson)	D10 DL10 DF10ST DFL10ST
Type Microliter pipette (design)	Air cushion pipette
Manufacturer	Gilson
Measurement deviation systematic in %	1 %
Measurement deviation systematic in µl	0.2 µl
Measurement deviation random in %	0.3 %
Measurement deviation random in µl	0.06 µl
Type pipette	Microliter pipettes
Type certificate	Declaration of Conformity (standardization) IVD Directive 98/79/EC

