



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson PIPETMAN® L P200L microliter pipette, metal ejector, volume 20–200 µL

€327.00
plus VAT &
Shipping

Product Images



Description

Single-channel pipette, variable type

More productivity more convenience

- ergonomic design
- light in weight
- with metal ejector
- 3 tip ejector positions - ideal for left- and right-handed users

GLP Features

- Volume lock/unlock mechanism
- unique 2D identification code
- traceable serial number
- color coded push buttons

Additional Information

No.	CF-0628
Manufacturer (Brand)	Gilson
old neoLab article no.	7-3404
EAN	4058072223007
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Volume MAX	200 µl
Volume MIN	20 µl
Operation via	Knobs/Buttons
Microliter pipettes features	with certification with GLP functions with color code ergonomically shaped IVD compliant
Number of channels	1
Designation tip (Gilson)	D10 DL10 DF10ST DFL10ST
Energy supply	without power supply
Color code	yellow
Functionality	mechanical
Manufacturer	Gilson
Measurement deviation systematic in %	0.8 %
Measurement deviation systematic in µl	1.6 µl
Measurement deviation random in %	0.15 %
Measurement deviation random in µl	0.3 µl
Model Microliter pipette/ Dispenser	Pipetman L P200L
Division	0.2 µl
Type Microliter pipette (design)	Air cushion pipette
Type pipette	Microliter pipettes
Type certificate	Declaration of Conformity (standardization) IVD Directive 98/79/EC

